

# Test Report



No: 220055202

Verification Website: [www.gtgc.net.cn](http://www.gtgc.net.cn)

Verification Code: BYRE-7483-54

Issue Date: 2022-03-23

Scan to download



Applicant: Hubei MAIYUE non woven products Co., LTD  
Address: Factory building no. 12-1, Zhanpeng Road, Zhongling Industrial Park, PENGCHANG  
Town, Xiantao City

Information confirmed by applicant:

Disposable face masks (KF94)

Quantity: 10 pieces

Basis of judgement:

GB 2626-2019 <Respiratory protection-Non-powered air-powered air-purifying particle respirator>

Date Received/Date Test Started: 2022-03-14

Conclusion:

Inspiratory resistance

M

Expiratory resistance

M

Note: "M"-Meet the standard's requirement "F"-Fail to meet the standard's requirement  
"---"-No comment

Remark:

This report is the english translation version of the report 220055201.

All the tested items are tested under the standard condition (except for indication).

Copies of the report are valid only re-stamped.

The experiment was carried out at No. 1, Zhujiang Road, Panyu District, Guangzhou, Guangdong, P. R. China.

Approved By:


JuXiang Zhong Engineer

*Ju Xiang Zhong*



# Test Report

No: 220055202



# Test Report

No: 220055202

---

## Inspiratory resistance (Pa)

M

(GB 2626-2019 6.5, Head size: Medium)

Untreated sample:

1# 140.1

2# 149.5

Pretreatment sample:

1# 136.2

2# 138.7

Requirement

≤210

(GB 2626-2019)

---

## Expiratory resistance (Pa)

M

(GB 2626-2019 6.6, Head size: Medium)

Untreated sample:

1# 132.6

2# 130.7

Pretreatment sample:

1# 127.8

2# 125.3

Requirement

≤210

(GB 2626-2019)

—End of Report—

Page

3/3



#### Note

1. The report is invalid without authorized stamp.
2. Copies of this report are invalid without authorized stamp.
3. Any dispute should be raised within 15 days after receiving the report.
4. The result is only valid for the tested sample.
5. The results of unapproved items are for reference only.
6. This report is invalid if altered.
7. Our company does not accept any responsibility for the authenticity of the information supplied by customers (including sample information).
8. The report shall not be duplicated separately or partly, without prior written permission approved by GTTC, except duplicated in full version.